

# HemoCue Safety Lancet

In order to ensure precise test results from capillary blood samples it is vital to have a spontaneous flow of blood. The HemoCue Safety Lancet makes a sufficiently deep incision of 2.25 mm. The device incorporates a particularly sharp needle and split-second incision to reduce the incidence of pain. To ensure that the lancet cannot be re-used, the needle is completely concealed inside the plastic casing after use.

## Data

Material: 23 gauge surgical steel needle

Depth: 2.25 mm

Package: 200 lancets in a cardboard box

